biber

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Read the instructions for use completely and carefully!

Technical data


- Material: polypropylene
- Weight: 205 grams
- Size: 20 cm x 23 cm x 8 cm

Not suitable for children
below 3 years

Components

1. Main saddle
2. Flaps
3. Steel spring

Introduction

The muscles of the pelvic floor are those muscle strands which run from the pubic bone to the coccyx and form a kind of loop around the genital area: around the urethra, vagina and anus in women or the anus and penile root in men – musculus pubococcygeus (lat. for pubococcygeal muscle) – PC for short.

- For women especially, targeted training of these muscles is important for a variety of reasons.
- A trained pelvic floor enhances pleasurable sensations and the ability to experience orgasm during sex. The better you are able to tense and relax these muscles, the more passionate will be your experience.
- Women can take action themselves to prevent prolapse of the uterus. Pelvic floor gymnastics are an ideal preparation for (natural) birth. Pregnancy and birth, in particular, can lead to overstretching of the pelvic floor muscles. Targeted training can return the muscles which have been intensely stretched during birth to their former shape. This will also bring back enjoyment of sex following a birth.
- Targeted pelvic floor training also helps with muscle weakening, especially during menopause.
- Trained pelvic floor muscles can prevent or reduce unwanted urinary incontinence in both women and men. Men are especially likely to suffer urinary incontinence following prostate operations.

It is never too late to start your pelvic floor training. Whether you start at 20 or at 70, what matters most is that you do the exercises regularly. Make gymnastics part of your daily routine.

Instructions for safe use of the product

- Keep these instructions for later!
- Immediately following operations, especially in the pelvic floor area, you should only start the training after consulting your treating physician. Your physician can decide how long following an operation you should go easy on your muscles.
- In order to prevent trapping skin, only use the pelvic floor trainer above clothing. Side effects that affect the entire organism are not yet known.
- People with limited physical, sensory or mental abilities should only use the pelvic floor exercise under supervision.
- If your symptoms worsen or you notice unusual changes, consult your doctor immediately.
- For additional questions consult your doctor!
- Due to the material composition, it can be worn directly on the skin as long as there are no external or internal injuries.
- Examine the appliance frequently. If there is any damage visible. If the supply cord is damaged, it must be replaced by the manufacturer or a qualified person, to avoid a hazard.
- All serious incidents involving the product must be reported to the manufacturer and the responsible authorities of the Member State where the user and/or patient is located.

Contraindications / Not suitable

Allergic reactions to the material

For Germany
Federal Institute for Drugs and Medical Devices
Medical Devices Division
Kurt-Georg-Kiesinger-Allee 3
53175 Bonn
<https://www2.bfarm.de/medprod/mpsv/>

For Austria
Federal Ministry of Health
Dept. Pharmaceuticals and Medical Devices
Radetzkystrasse 2, 1030 Vienna
Email: meddev@bmgf.gv.at

For Switzerland
Swissmedic, Swiss Agency for
Therapeutic Products
Medical Devices Division
Hallerstrasse 7, 3012 Bern, Switzerland

Intended use


The pelvic floor trainer is used to strengthen pelvic floor muscles. The aim is to train your muscles a little more every day without exerting too much of a strain or causing pain.

Applicable to

- Disc discomfort
- Prevention of uterine subsidence
- Reduces or prevents urinary incontinence

For private use.

Use


- Wear comfortable clothing and do some basic exercises to warm up before starting your training with the pelvic floor trainer.
- Place the trainer between your thighs, several centimetres below your pelvic floor. Do NOT place it right up against your pelvic floor.
- Using the flaps, centre the pelvic floor trainer on the front or back of your thighs, whichever you find the more comfortable.
- Make sure you keep the right posture. Stand up straight with your feet close together and your knees slightly bent.
- Try to keep your knees, hips and shoulders in a line and press the flaps together using your thighs, buttocks and pelvic floor muscles.
- Beginners should press and hold the flaps for 3 seconds and then rest for 7 seconds. Try to manage six reps.

You will initially find it difficult to close the flaps. That's fine because the exercise will help you activate your major muscles even if you only close the flaps halfway.

- Increase use and the hold time gradually. Start by using the trainer three times a week and build this up to twice a day. Give your muscles enough time to recover. Avoid over-exerting your muscles.
- Trained users should manage 10 reps. Tense for 25 seconds, relax for 5 seconds.
- If you train regularly with the right positioning, your muscle strength will increase automatically.

Make sure to always push the flaps together as far as possible. The training is efficient and the results achieved can vary from person to person.

Other applications includes training your inner thigh muscles or chest muscles:

Cleaning and storage

- Store the device in a clean, dry place.
- Do not expose the pelvic floor trainer to direct sunlight for a long time.
- Clean the trainer after each use by wiping it with a dry or moist cloth.

Warranty

- The warranty is valid for 24 months following the purchase date and covers manufacturer's faults relating to material and quality. Liability for damage cannot be accepted if the device is used for a purpose other than intended or if it is used incorrectly.
- This warranty is only valid if the instructions in the manual are observed. If damage is caused by force, misuse or inappropriate use or servicing by unauthorised parties, the warranty ceases to apply.

Disposal

This product can be safely disposed of in household waste.

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Biber Umweltprodukte Versand GmbH
Hydas GmbH & Co. KG, Hirzenhainer Straße 3,